

MX-PRO® R3

x-frame ambulance cot

stryker®



EMS Equipment

3800 E. Centre Ave
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 784 4336

www.ems.stryker.com

Standard Features

- Transport at load height capability
- Positive action height adjustment
- Easy to use release handle design
- Color-coded controls
- Dampened action during hot drops
- High visibility powder-coated frame
- Light-weight, rugged aluminum construction
- Scientifically-optimized lift bar design
- Lower lifting bar
- Seven height positions
- Lift-capable safety bar
- Integrated bumper system
- Perforated litter surface
- One-hand release breakaway head section with safety bar
- Floor-mounted safety hook
- One-hand release fold-down siderails
- One-hand release, infinite positioning, lift assist pneumatic backrest
- Oversized wheels with sealed caster and wheel bearings
- Reflective labeling
- Sealed bolster mattress
- Shock position
- Two lap belts and one four-point shoulder restraint
- Single wheel lock

Optional Features

- Heavy duty two- or three-stage IV poles (patient right or left)
- Base lift bar
- Hard or soft base storage
- High height kit for 33"
- Defibrillator platform
- Pocketed or non-pocketed head-end storage
- Height limit kit*
- Permanent or removable oxygen bottle holders (head or foot end)
- Pull handle
- Sealed flat mattress
- Premium mattress
- Dual wheel lock
- X-frame guards
- Head extension
- Pillow
- Rigid head-end storage tray
- Equipment hook
- Backrest storage pouch

Specifications

• Model Number	6082	
• Height¹		
Position 1	13"	(33 cm)
Position 2	20.5"	(52 cm)
Position 3	25"	(64 cm)
Position 4	29"	(74 cm)
Position 5	32"	(81 cm)
Position 6	35"	(89 cm)
Position 7	38"	(97 cm)
• Length		
Standard	80"	(203 cm)
Minimum	62"	(157 cm)
• Width	23"	(58 cm)
• Weight²	81 lbs	(37 kg)
• Wheels		
Diameter	6"	(15 cm)
Width	2"	(5 cm)
• Articulation		
Backrest	0° - 73°	
Shock Position	+14°	
• Maximum Weight Capacity	650 lbs	(295 kg)
• Minimum Operators Required		
Occupied Cot	2	
Unoccupied Cot	1	
• Recommended Fastener System		
Floor Mount	Model 6370 or 6377	
Wall Mount	Model 6371	
• Recommended Loading Height*	Up To 33"	(84 cm)

Warranty

- One year parts, labor and travel or two years parts

*Height limit kit recommended for load heights less than 28" (71 cm).

High height kit recommended for load heights more than 32" (81 cm).

¹Height measured from bottom of mattress at seat section to ground level.

²Cot is weighed without mattress and restraints.

Stryker reserves the right to change specifications without notice.

In-service video included with every order.

The MX-PRO® R3 is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822D.

The MX-PRO® R3 is designed to be compatible with competitive cot fastener systems.