

Power-PRO™ XT

powered ambulance cot

stryker®



EMS Equipment

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100 f: 866 795 2233
toll free: 800 784 4336

www.ems.stryker.com

Standard Features

- Automatic in-cot fastener shut-off
- Manual back-up system
- Automatic high-speed retract
- Battery-powered hydraulic lift system
- Preventative maintenance hour meter
- Battery power indicator
- Settable load height with jog function
- Color-coded controls
- High visibility powder-coated frame
- Lightweight, rugged aluminum construction
- Scientifically optimized lift bar and control design
- Lower lifting bar
- Lift-capable safety bar
- Integrated bumper system
- Retractable head section
- Floor-mounted safety hook
- One-hand release, fold down side rails
- One-hand release, infinite positioning, lift assist pneumatic backrest
- Oversized wheels with sealed caster and wheel bearings
- Sealed bolster mattress
- Shock or flat leg positioning
- Two lap belts and one four-point shoulder restraint
- Single wheel lock
- Pull handle
- X-frame guards
- Power washable

Optional Features

- Heavy duty two- or three-stage IV poles (patient right or left)
- Permanent or removable oxygen bottle holders (head or foot)
- Sealed flat mattress
- Dual wheel locks
- Head extension
- Pillow
- Equipment hook
- Backrest storage pouch
- Head end storage flat
- Defibrillator platform
- Transfer flat
- Base storage net
- Knee gatch
- SMRT Power System (12V DC, 120V AC and 240V AC options available)
- SMRT Charger mounting bracket
- DEWALT® Battery System (12V, 110V and 220V models available)

Specifications

Model Number	6500
Height¹ (infinite height positioning between lowest and highest position)	
Lowest Position	14" (36 cm)
Highest Position	41.5" (105 cm)
Length	
Standard	81" (206 cm)
Minimum	63" (160 cm)
Width	23" (58 cm)
Weight²	125 lbs (57 kg)
Wheels	
Diameter	6" (15 cm)
Width	2" (5 cm)
Articulation	
Backrest	0 – 73°
Shock Position	+15°
Knee Gatch	30°
Maximum Weight Capacity³	700 lbs (318 kg)
Minimum Operator Required	
Occupied Cot	2
Unoccupied Cot	1
Recommended Fastener System	
Floor Mount	Model 6370 or 6377
Wall Mount	Model 6371
Recommended Loading Height⁴	Up to 36" (91 cm)

Warranty

- Two-year parts, labor and travel
- One-year soft goods
- Three-year X-frame components
- Three-year limited powertrain
- Lifetime on all welds*

Extended warranties available

¹ Height measured from bottom of mattress, at seat section, to ground level.

² Cot is weighed with 1 battery pack, without mattress and restraints.

³ 700 lb weight capacity with an unassisted lift capacity of 500 lbs.

⁴ Load wheel height can be set between 26" and 36".

Stryker reserves the right to change specifications without notice.

* 7-year service life.

In-service video included with every order.

The Power-PRO is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822E.

The Power-PRO is designed to be compatible with competitive cot fastener systems.

DEWALT is a registered trademark of DEWALT Industrial Tool Company.

Patents pending.



Copyright © 2009 Stryker
PC/RRP 5M 1009, 08-SEM-2328
Mkt Lit-152 Rev C